

Scaler marked with Color Coding Rings

(IMS-1280L, IMS-1286L, IMS-1281, IMS-1287, IMS-12810, IMS-1287L, IMS-12810L, IMS-1288, IMS-12811, IMS-1288L, IMS-12811L, IMS-1289, IMS-1281L, IMS-1289L, IMS-1282, IMS-1282L, IMS-1283, IMS-1283L, IMS-1284, IMS-1284L, IMS-1234, IMS-1285, IMS-1285L, IMS-1280, IMS-1286)

Processing:

- For reprocessing, attached **Color Coding Rings** do not have to be removed
- Validation reports:
 - Project Number: 00418-1
Examination of an Automated Cleaning Process using quantitative Detection of Protein and Hemoglobin and the Radionuclide Method
07-Jun-2018
 - Project Number: 00418-2
Examination of a Manual Cleaning Process using quantitative Detection of Protein and Hemoglobin and the Radionuclide Method
07-Jun-2018
 - Project Number: 25517-1
Validation of a Sterilization Process using Steam Sterilization in Pre-vacuum Mode
Method MD 4.0: Sterilization validation of medical products with moist heat
12-Dec-2018
 - Project Number: 25517-2
Validation of a Sterilization Process using Steam Sterilization in Pre-vacuum Mode
Method MD 4.0: Sterilization validation of medical products with moist heat
12-Dec-2018